


Product Name: Recombinant Human UPB1 Protein
Catalog #: PEH7986



Summary

Name	BUP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UBR1
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAGAEWKSLE ECLEKHLPLP DLQEVKRVLY GKELRKLDLP REAFEASRE DFELQGYAFE AAEEQLRRPR IVHVGLVQNR IPLPANAPVA EQVSALHRR KAIVEVAAMC GVNII CFQEA WTMPFAFCTR EKLPWTEFAE SAEDGPTTRF CQKLAKNHDM VVVSPILERD SEHGDVLWNT AVVISNSGAV LGKTRKNHIP RVGDFNESTY YMEGNLGHPV FQTQFGRIAV NICYGRHHPL NWLMYSINGA EIIFNPSATI GALSESLWPI EARNAAIANH CFTCAINRVG TEHFPNEFTS GDGKKAHQDF GYFYGSSYVA APDSSRTPGL SRSRDGLLVA KLDLNLCCQQV NDVWNFKMTG RYEMYARELA EAVKSNYSPT IVKE

image



Background

Alternative Names	Beta-ureidopropionase (EC 3.5.1.6) (BUP-1) (Beta-alanine synthase) (N-
--------------------------	--

Product Name: Recombinant Human UPB1 Protein
Catalog #: PEH7986



carbamoyl-beta-alanine amidohydrolase)

Note

For Research Use Only , Not for Diagnostic Use