


Product Name: Recombinant Human UBA2 Protein
Catalog #: PEH7987



Summary

Name	SAE2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UBT2
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PGCTIRNTPS EPIHCIVWAK YLFNQLFGEE DADQEVSPDR ADPEAAWEPT EAEARARASN EDGDIKRIST KEWAKSTGYD PVKLF TKLFK DDIRYLLTMD KLWRKRKPPV PLDWAEVQSQ GEETNASDQQ NEPQLGLKDQ QVLDVKSYPAR LFSKIETLR VH LAEKGDGA ELIWDKDDPS AMDFV TSAAN LRMHIFSMNM KSRFDIKSMA

image



Background

Alternative Names	SUMO-activating enzyme subunit 2 (EC 2.3.2.-) (Anthracycline-associated resistance ARX) (Ubiquitin-like 1-activating enzyme E1B) (Ubiquitin-like modifier-activating enzyme 2)
--------------------------	--

Product Name: Recombinant Human UBA2 Protein
Catalog #: PEH7987



Note

For Research Use Only , Not for Diagnostic Use