


Product Name: Recombinant Human NXF1 Protein
Catalog #: PEH7988



Summary

Name	NXF1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UBU9
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GRGGSGIRSS RLEEDGDVA MSDAQDPRV RYNPYTTRPN RRGDTWHDRD RIHVTVRRDR APPERGGAGT

image



Background

Alternative Names Nuclear RNA export factor 1 (Tip-associated protein) (Tip-associating protein) (mRNA export factor TAP)

Note

For Research Use Only , Not for Diagnostic Use