


Product Name: Recombinant Human CNOT8 Protein
Catalog #: PEH7991



Summary

Name	CNOT8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UFF9
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LQFQKHEEEG IDTLHFAELL MTSGVVLCDN VKWLSFHSGY DFGYMKLLT DSRLPEEEHE FFHILNLFFP SIYDVKYLMK SCKNLKGGLQ EVADQLDLQR IGRQHQAGSD SLLTGMAFFR MKELFFEDSI DDAKYCGRLY GLGTGVAQKQ NEDVDSAQEK MSILAIINNM

image



Background

Alternative Names CCR4-NOT transcription complex subunit 8 (EC 3.1.13.4) (CAF1-like protein) (CALIFp) (CAF2) (CCR4-associated factor 8) (Caf1b)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CNOT8 Protein
Catalog #: PEH7991

