


Summary

| | |
|------------------------------|---|
| Name | TCF20_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q9UGU0 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 58 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | ISIPSKRQS QDVKNSSTED KGRLLHSSKE GADKAFNSYA HLSHSQDIKS IPKRDSSKDL PPSDRNCPA VTLTSPAKTK ILPPRKGRGL KLEAIVQKIT SPNIRRSASS NSAEAGGDTV TLDDILSLKS GPPEGGSVAV QDADIEKRKG EVASDLVSPA NQELHVEKPL PRSSEEWGRS VDDKVKTETH AETVTAGKEP PGAMTSTTSQ KPGSNQGRPD GSLGGTAPLI FPDSKNVPPV GILAPEANPK AEEKENDTVT ISPKQEGFPP KGYFPSGKKK GRPIGSVVKQ KKQQQPPPPP PQPPQIPEGS ADGEPKPKKQ RQRRERRKPG AQPRKRKTKQ AVPIVEPQEP EIKLKYATQP LDKTDAKNKS FYPYIHVVNK CELGAVCTII NAEETEQTKL VRGRKGQRSL TPPPSSTESK ALPASSFMLQ GPVVTESSVM GHLVCCLCGK WASYRNMGDL FGPFPQDYA ATLPKNPPPK RATEMQSKVK VRHKSASNGS KTDTEEEEEEQ QQQQKEQRSL |

image



Background

Product Name: Recombinant Human TCF20 Protein
Catalog #: PEH7992



Alternative Names

Transcription factor 20 (TCF-20) (Nuclear factor SPBP) (Protein AR1)
(Stromelysin-1 PDGF-responsive element-binding protein) (SPRE-binding
protein)

Note

For Research Use Only , Not for Diagnostic Use