


Product Name: Recombinant Human TBK1 Protein
Catalog #: PEH7994



Summary

Name	TBK1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q9UHD2
Host	E.coli
Species	Human
Molecular Mass	53 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EIYYQFKKDK AERRLAYNEE QIHKFDKQKL YYHATKAMTH FTDECVKKYE AFLNKSEEWI RKMLHLRQQL

image



Background

Alternative Names Serine/threonine-protein kinase TBK1 (EC 2.7.11.1) (NF-kappa-B-activating kinase) (T2K) (TANK-binding kinase 1)

Note

For Research Use Only , Not for Diagnostic Use