


Product Name: Recombinant Human USP25 Protein
Catalog #: PEH7995



Summary

Name	UBP25_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UHP3
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LNQLREITGI NDTQILQQAL KDSNGNLELA VAFLTAKNAK TPQQEETYY QTALPGNDRY ISVGSQADTN VIDLTGDDKD DLQRAIALSL AESNRAFRET GITDEEQAIS RVLEASIAEN KACLRTPTE VWRDSRNPYD

image



Background

Alternative Names	Ubiquitin carboxyl-terminal hydrolase 25 (EC 3.4.19.12) (Deubiquitinating enzyme 25) (USP on chromosome 21) (Ubiquitin thioesterase 25) (Ubiquitin-specific-processing protease 25)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human USP25 Protein
Catalog #: PEH7995

