


Product Name: Recombinant Human XPO7 Protein
Catalog #: PEH7997



Summary

Name	XPO7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q9UIA9
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VHILRDGLE DPLEDTGLVQ QQLDQLSTIG RCEYEKTCAL LVQLFDQSAQ SYQELLQSAS ASPMDIAVQE GRLTWLVYII GAVIGGRVSF ASTDEQDAMD GELVCRVLQL MNLTDSRLAQ AGNEKLELAM LSFFEQFRKI YIGDQVQKSS KLYRRLSEVL GLNDETMVLS VFIGKIITNL KYWGRCEPIT SKTLQLLNDL

image



Background

Alternative Names Exportin-7 (Exp7) (Ran-binding protein 16)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human XPO7 Protein
Catalog #: PEH7997

