


Product Name: Recombinant Human DCTN4 Protein
Catalog #: PEH8004



Summary

Name	DCTN4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UJW0
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QKEIKIEPAQ AVDEVEPLPE DYYTRPVNLT EVTTLQQRLL QPDFQPVCAS QLYPRHKHLL IKRSLRCRKC EHNLSKPEFN PTSIKFKIQL VAVNYIPEVR IMSIPNLRYM KESQVLLTLT NPVENLTHVT LFECEEGDPD DINSTAKVVV PPKELVLAGK DAAAEYDELA EPQDFQDDPD IIAFRKANKV GIFIKVTPQR EEGEVTVCFK MKHDFKNLAA PIRPIEESDQ

image



Background

Alternative Names Dynactin subunit 4 (Dyn4) (Dynactin subunit p62)

Note

Product Name: Recombinant Human DCTN4 Protein
Catalog #: PEH8004



For Research Use Only , Not for Diagnostic Use