


Product Name: Recombinant Human PTBP2 Protein
Catalog #: PEH8006



Summary

Name	PTBP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UKA9
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GVSAGGNTVL LVSNLNEEMV TPQSLFTLFG VYGDVQRVKI LYNKKDSALI QMADGNQSQL AMNHLNGQKM YGKIIRVTLK KHQTVQLPRE GLDDQGLTKD FGNSPLHRFK KPGSKNFQNI FPPSATLHLS NIPPSVAEED LRTL FANTGG TVKAFKFFQD HKMALLQMAT VEEAIQALID LHNYNLGENH HLRVSFSKST

image



Background

Alternative Names Polypyrimidine tract-binding protein 2 (Neural polypyrimidine tract-binding protein) (Neurally-enriched homolog of PTB) (PTB-like protein)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PTBP2 Protein
Catalog #: PEH8006

