



Product Name: Recombinant Human IKZF3 Protein
Catalog #: PEH8008



Summary

Name	IKZF3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UKT9
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KSTQEQSVPA ESAAVLNDYS LTKSHEMENV DSGEGPANED EDIGDDSMKV KDEYSERDEN VLKSEPMGNA EEPEIPYSYS REYNEYENIK LERHVVSFDS SRPTSGKMNC DVCGLSCISF NVLMVHKRSH TGERPFQCNQ CGASFTQKGN LLRHIKLHTG EKPFKCHLCN YACQRRDALT GHLRTHSVEK PYKCEFCGRS YKQRSSLEEH KERCRTFLQS TDPGDTASAE ARHIKAEMGS ERALVLDRLA SNVAKRKSSM

image



Background

Alternative Names Zinc finger protein Aiolos (Ikaros family zinc finger protein 3)

Product Name: Recombinant Human IKZF3 Protein
Catalog #: PEH8008



Note

For Research Use Only , Not for Diagnostic Use