


Product Name: Recombinant Human ASAP1 Protein
Catalog #: PEH8013



Summary

Name	ASAP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9ULH1
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PTVDIVNQAG ETALDIAKRL KATQCEDLLS QAKSGKFNPH VHVEYEWNL QEEIDESDDD LDDKPSPKK ERSRPQSFC HSSISPQDK LALPGFSTPR DKQRLSYGAF

image



Background

Alternative Names	Arf-GAP with SH3 domain, ANK repeat and PH domain-containing protein 1 (130 kDa phosphatidylinositol 4,5-bisphosphate-dependent ARF1 GTPase-activating protein) (ADP-ribosylation factor-directed GTPase-activating protein 1) (ARF GTPase-activating protein 1) (Development and differentiation-enhancing factor 1) (DEF-1) (Differentiation-enhancing factor 1) (PIP2-dependent ARF1 GAP)
--------------------------	--

Product Name: Recombinant Human ASAP1 Protein
Catalog #: PEH8013



Note

For Research Use Only , Not for Diagnostic Use