


Product Name: Recombinant Human VPS4A Protein
Catalog #: PEH8019



Summary

Name	VPS4A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UN37
Host	E.coli
Species	Human
Molecular Mass	50 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MTTSTLQKAI DLVTKATEED KAKNYEEALR LYQHAVEYFL HAIKYEAHSD KAKESIRAKC VQYLDRAEKL KDYLRSKEKH GKPPVKENQS EGKGSDDSE GDNPEKKKLQ EQLMGAVVME KPNIRWNDVA GLEGAKEALK EAVILPIKFP HLFTGKRTPW RGILLFGPPG TGKSYLAKAV ATEANNSTFF SVSSDLMSK WLGESEKLVK NLFELARQHK PSIIFIDEVD SLCGSRNENE SEAARRIKTE FLVQMQGVGN NNDGTLVLGA TNIPWVLDSA IRRFEKRIY IPLPEEAARA QMFRLHLGST PHNLTDANIH ELARKTEGYS GADISIIVRD SLMQPVRKVQ SATHFKKVCG PSRTNPMMI DLLTPCSPG DPGAMEMTWM DVPGDKLLP VVCMSDMLRS LATRPTVNA DLLKVKKFS EDFGQES

image



Background

Product Name: Recombinant Human VPS4A Protein
Catalog #: PEH8019



Alternative Names

Vacuolar protein sorting-associated protein 4A (EC 3.6.4.6) (Protein SKD2)
(VPS4-1) (hVPS4)

Note

For Research Use Only , Not for Diagnostic Use