


Product Name: Recombinant Human SNX6 Protein
Catalog #: PEH8021



Summary

Name	SNX6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UNH7
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MVKSADGVIV SGVKDVDDFF EHRTFLEY HNRVKDASAK SDRMTRSHKS AADDYNRIGS SLYALGTQDS TDICKFFLKV SELFDKTRKI EARVSADEDL KLSDLLKYYL RESQAAKDLL YRRSRSLVDY ENANKALDKA RAKNKDVLQA ETSQQLCCQK FEKISESAKQ ELIDFKTRRV AAFRKNLVEL AELELKHAKG NLQLLQNCLA

image



Background

Alternative Names Sorting nexin-6 (TRAF4-associated factor 2) [Cleaved into: Sorting nexin-6, N-terminally processed]

Note

Product Name: Recombinant Human SNX6 Protein
Catalog #: PEH8021



For Research Use Only , Not for Diagnostic Use