


Product Name: Recombinant Human TIMELESS Protein
Catalog #: PEH8023



Summary

Name	TIM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UNS1
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DQAIVSGNVP SSPEEVEAVW PALAEQLQCC AQNSELSMDS VVPFDAASEV PVVEEQRAEAM VRIQDCLLAG QAPQALTLR SAREVWPEGD VFGSQDISPE EEIQLLKQIL SAPLPRQQGP EERGAE EEEEE EEEEEEEELQ

image



Background

Alternative Names Protein timeless homolog (hTIM)

Note

For Research Use Only , Not for Diagnostic Use