


**Product Name: Recombinant Human BSN Protein**  
**Catalog #: PEH8025**



## Summary

<b>Name</b>	BSN_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9UPA5
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	53 kDa
<b>Buffer</b>	Lyophilized from Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	ELEMESLTGS PEDRSRGEHS STLPASTPSY TSGTSPTSLS SLEEDSDSSP SRRQRLEEAK QQRKARHRSH GPLLPTIEDS SEEEELREEE ELLREQEKMR EVEQQRIRST ARKTRRDKEE LRAQRRRERS KTPPSNLSPI EDASPTTELR QAAEMEELHR SSCSEYSPSP SLDSEAEALD GGPSRLYKSG SEYNLPTFMS LYSPTETPSG SSTTPSSGRP LKSAEEAYEE MMRKAELLQR QQGQAAGARG PHGGPSQPTG PRGLGSFEYQ DTTDREYGQA AQPAAEGTPA SLGAAVYEEI LQTSQSIVRM RQASSRDLAF AEDKKKEKQF LNAESAYMDP MKQNGGPLTP GTSPTQLAAP VSFSTPTSSD SSGGRVIPDV RVTQHFAKET QDPLKLHSSP ASPSSASKEI GMPFSQPGT PATTAVAPCP AGLPRGYMTP ASPAGSERSP SPSSTAHSYG HSPTTANYGS

## image



## Background

**Product Name: Recombinant Human BSN Protein**  
**Catalog #: PEH8025**



---

**Alternative Names**                      Protein bassoon (Zinc finger protein 231)

**Note**

For Research Use Only , Not for Diagnostic Use