

Product Name: Recombinant Human SATB2 Protein
Catalog #: PEH8027



Summary

Name	SATB2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UPW6
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IQQEMKRAKV SQALFAKVAA NKSQGWLCLE LRWKENPSPE NRTLWENLCT IRRFLNLPQH ERDVIYEEES RHHHSERMQH VVQLPPEPVQ VLHRQQSQPA KESSPPREEA PPPPPPTEDS CAKKPRSRTK ISLEALGILQ SFIHDVGLYP DQEIHTLSA QLDLPKHTII KFFQNQRYHV KHHGKLKEHL GSAVDVAEYK DEELLTESEE NDSEEGSEEM YKVEAEEENA

image



Background

Alternative Names DNA-binding protein SATB2 (Special AT-rich sequence-binding protein 2)

Note

Product Name: Recombinant Human SATB2 Protein
Catalog #: PEH8027



For Research Use Only , Not for Diagnostic Use