


## Summary

|                              |  |
|------------------------------|--|
| <b>Name</b>                  | DYN3_HUMAN   |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE  |
| <b>Endotoxin level</b>       | Please contact with the lab for this information   |
| <b>Tag</b>                   | C-His  |
| <b>Accession #</b>           | Q9UQ16   |
| <b>Host</b>                  | E.coli   |
| <b>Species</b>               | Human  |
| <b>Molecular Mass</b>        | 52 kDa   |
| <b>Buffer</b>                | Lyophilized from Phosphate buffered saline   |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.  |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.   |
| <b>Sequence</b>              | LPQIAVGGGQ SAGKSSVLEN FVGRDFLPRG SGIVTRRPLV LQLVTSKAEY<br>AEFLHCKGKK FTFDFEVRLE IEAETDRVTG MNKGISSIPI NLRVYSPHVL<br>NLTLIDLPGI TKVPVGDQPP DIEYQIREMI MQFITRENCL ILAVTPANTD<br>LANSALKLA KEVDPQGLRT IGVITKLDLM DEGTDARDVL ENKLLPLRRG<br>YGVVNRSQK DIDGKKDIKA AMLAERKFFL SHPAYRHIAD RMGTPHLQKV<br>LNQQLTNHIR DTLPNFRNKL QGQLLSIEHE VEAYKNFKPE DPTRKTKALL<br>QMVQQFAVDF EKRIEGSGDQ VDTLELSGGA KINRIFHERF PFEIVKMEFN<br>EKELRREISY AIKNIHGIRT GLFTPDMAFE AIVKKQIVKL KGPSLKSVDL<br>VIQELINTVK KCTKLANFP RLCEETERIV ANHIREREGK |

## image



## Background

**Product Name: Recombinant Human DNMT3 Protein**  
**Catalog #: PEH8028**



---

**Alternative Names**

Dynamin-3 (EC 3.6.5.5) (Dynamin, testicular) (T-dynamin)

**Note**

For Research Use Only , Not for Diagnostic Use