


Product Name: Recombinant Human CRYL1 Protein
Catalog #: PEH8033



Summary

Name	CRYL1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y2S2
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IVGSGVIGRS WAMLFASGGF QVKLYDIEQQ QIRNALENIR KEMKLLEQAG SLKGSLSVEE QLSLISGCPN IQEAVEGAMH IQECVPEDLE LKKKIFAQLD SIIDDRVILS SSTSCLMPK LFAGLVHVKQ CIVAHPVNPP YYIPLVELVP HPETAPTTVD RTHALMKKIG QCPMRVQKEV AGFVLNRLQY AIISEAWRLV EEGIVSPSDL DLVMSEGLGM RYAFIGPLET MHLNAEGMLS YCDRYSEGIK HVLQTFGPIP EFSRATAEKV NQDMCMKVPD DPEHLAARRQ WRDECLMRLA

image



Background

Alternative Names	Lambda-crystallin homolog (EC 1.1.1.45) (L-gulonate 3-dehydrogenase) (Gul3DH)
--------------------------	---

Product Name: Recombinant Human CRYL1 Protein
Catalog #: PEH8033



Note

For Research Use Only , Not for Diagnostic Use