


Product Name: Recombinant Human UBE2J1 Protein
Catalog #: PEH8036



Summary

Name	UB2J1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y385
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AKKSQDFCCE GCGSAMKDVL LPLKSGSDSS QADQEAKELA RQISFKAENV SSGKTISESD LNHSFSLTDL QDDIPTTFQG ATASTSYGLQ NSSAASFHQP TQPVAKNTSM SPRQRRRAQQQ SQRRLSTSPD VIQGHQPRDN

image



Background

Alternative Names	Ubiquitin-conjugating enzyme E2 J1 (EC 2.3.2.23) (E2 ubiquitin-conjugating enzyme J1) (Non-canonical ubiquitin-conjugating enzyme 1) (NCUBE-1) (Yeast ubiquitin-conjugating enzyme UBC6 homolog E) (HsUBC6e)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human UBE2J1 Protein
Catalog #: PEH8036

