


Product Name: Recombinant Human REXO2 Protein
Catalog #: PEH8038



Summary

Name	ORN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y3B8
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VREGGAAMAA GESMAQRMVW VDLEMTGLDI EKDQIEMAC LITDSDLNIL AEGPNLIKQ PDELLDSMSD WCKEHHGKSG LTKAVKESTI TLQQAIEYFL SFVRQQTPPG LCPLAGNSVH EDKKFLDKYM PQFMKHLHYR IIDVSTVKEL CRRWYPEEYE FAPKKAASHR ALDDISESIK ELQFYRNNIF KKKIDEKKRK IIENGENEKT VS

image



Background

Alternative Names Oligoribonuclease, mitochondrial (EC 3.1.-.-) (RNA exonuclease 2 homolog)
(Small fragment nuclease)

Note

Product Name: Recombinant Human REXO2 Protein
Catalog #: PEH8038



For Research Use Only , Not for Diagnostic Use