


Product Name: Recombinant Human TSC22D4 Protein
Catalog #: PEH8041



Summary

Name	T22D4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y3Q8
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SGHLDSDDDS GSGSLVGIDN KIEQAMD LVK SHLMFAVREE VEVLKEQIRE LAERNAALEQ ENGLLRALAS

image



Background

Alternative Names	TSC22 domain family protein 4 (TSC22-related-inducible leucine zipper protein 2) (Tsc-22-like protein THG-1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use