


Product Name: Recombinant Human SAMHD1 Protein
Catalog #: PEH8043



Summary

Name	SAMH1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y3Z3
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LKADDYIEIT GAGGKKYRIS TAIDDMEAYT KLTDNIFLEI LYSTDPKPKD AREILKQIEY RNLFKYVGET QPTGQIKIKR EDYESLPKEV ASAKPKVLLD VKLKAEDFIV DVINMDYGMQ

image



Background

Alternative Names	Deoxynucleoside triphosphate triphosphohydrolase SAMHD1 (dNTPase) (EC 3.1.5.-) (Dendritic cell-derived IFNG-induced protein) (DCIP) (Monocyte protein 5) (MOP-5) (SAM domain and HD domain-containing protein 1) (hSAMHD1)
--------------------------	--

Note

Product Name: Recombinant Human SAMHD1 Protein
Catalog #: PEH8043



For Research Use Only , Not for Diagnostic Use