


Product Name: Recombinant Human GPSM3 Protein
Catalog #: PEH8045



Summary

Name	GPSM3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y4H4
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DGEQGPPQDE EGWPPPNSTT RPWRSAPPSP PPPGTRHTAL GPRSASLLSL QTELLDLVA EAQSRRL EEQ RATFYTPQNP SSLAPAPLRP LEDREQLYST ILSHQCQRME AQRSEPPLPP

image



Background

Alternative Names G-protein-signaling modulator 3 (Activator of G-protein signaling 4) (G18.1b) (Protein G18)

Note

For Research Use Only , Not for Diagnostic Use