


Product Name: Recombinant Human GMPPB Protein
Catalog #: PEH8046



Summary

Name	GMPPB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y5P6
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GTRLRPLTLS TPKPLVDFCN KPILLHQVEA LAAAGVDHVI LAVSYMSQVL EKEMKAQEQR LGIRISMSHE EEPLGTAGPL ALARDLLSET ADPFFVLNSD VICDFPFQAM VQFHRHHGQE GSILVTKVEE PSKYGVVVCE ADTGRIHRFV EKPQVFVSNK INAGMYILSP AVLQRIQLQP TSIEKEVPFI MAKEGQLYAM ELQGFWM DIG QPKDFLTGMC LFLQSLRQKQ

image



Background

Alternative Names	Mannose-1-phosphate guanylyltransferase beta (EC 2.7.7.13) (GDP-mannose pyrophosphorylase B) (GTP-mannose-1-phosphate guanylyltransferase beta)
--------------------------	---

Note

Product Name: Recombinant Human GMPPB Protein
Catalog #: PEH8046



For Research Use Only , Not for Diagnostic Use