


Product Name: Recombinant Human CDC42BPB Protein
Catalog #: PEH8047



Summary

Name	MRCKB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y5S2
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QRDLEHSLQM EAYERRIRRL EQEKLELSRK LQUESTQTVQS LHGSSRALSN SNRDKEIKKL NEEIERLKNK IADSNRLERQ LEDTVALRQE REDSTQRLRG LEKQHRVVRQ EKEELHKQLV EASERLKSQA KELKDAHQQR KLALQEFSEL NERMAELRAQ KQKVSRLRD KEEEMEVATQ KVDAMRQEMR RAEKLRKELE AQLDDAVAEA SKERKLREHS ENFCKQMESE LEALKVKQGG RGAGATLEHQ QEISKIKSEL EKKVLFYEEE LVRREASHVL EVKNVKKEVH DSESHQLALQ KEILMLKDKL ESKRERHNE MEEAVGTIKD

image



Background

Alternative Names	Serine/threonine-protein kinase MRCK beta (EC 2.7.11.1) (CDC42-binding protein kinase beta) (CDC42BP-beta) (DMPK-like beta) (Myotonic dystrophy
--------------------------	---



Product Name: Recombinant Human CDC42BPB Protein
Catalog #: PEH8047

kinase-related CDC42-binding kinase beta) (MRCK beta) (Myotonic
dystrophy protein kinase-like beta)

Note

For Research Use Only , Not for Diagnostic Use