


Product Name: Recombinant Human CLIC4 Protein
Catalog #: PEH8051



Summary

Name	CLIC4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y696
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RKPADLQNLA PGTHPPFITF NSEVKTDV NK IEEFLEE VLC PPKYLKLS PK HPESNTAGMD IFAKFSAYIK NSRPEANEAL ERLLKTLQK LDEYLNSPLP DEIDENSMED IKFSTRKFLD GNEMTLADCN LLPKLHIVKV VAKKYRNFDI PKEMTGIWRY LTNAYS RDEF TNTCPSDKEV EIAYS DVAKR

image



Background

Alternative Names Chloride intracellular channel protein 4 (Intracellular chloride ion channel protein p64H1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CLIC4 Protein
Catalog #: PEH8051

