


Product Name: Recombinant Human OAS3 Protein
Catalog #: PEH8053



Summary

Name	OAS3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q9Y6K5
Host	E.coli
Species	Human
Molecular Mass	56 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PNQKTPENSK SLNAVYPRAG SKPPSCPAPG PTGAASIVPS VPGMALDLSQ IPTKELDRFI QDHLKPSPQF QEQVKKAIDI ILRCLHENCV HKASRVSKGG SFGRGTDLRD GCDVELIIFL NCFTDYKDQG PRRAEILDEM RAQLESWWQD QVPSLSLQFP EQNVPEALQF QLVSTALKSW TDVSLPAPFD AVGQLSSGTK PNPQVYSRLL TSGCQEGEHK ACFAELRRNF MNIRPVKLKN LILLVKHWYR QVAAQNKGGK PAPASLPPAY ALELLTIFAW EQGCRQDCFN MAQGFRTVLG LVQQHQQLCV YWTVNYSTED PAMRMHLLGQ LRKPRPLVLD PADPTWNVGH GSWELLAQEA AALGMQACFL SRDGTSVQPW

image



Background

Alternative Names	2'-5'-oligoadenylate synthase 3 ((2-5')oligo(A) synthase 3) (2-5A synthase 3)
--------------------------	---

Product Name: Recombinant Human OAS3 Protein
Catalog #: PEH8053



(EC 2.7.7.84) (p100 OAS) (p100OAS)

Note

For Research Use Only , Not for Diagnostic Use