


Product Name: Recombinant Human USH1C Protein
Catalog #: PEH8054



Summary

Name	USH1C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y6N9
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RERLAEARQR ELQRQELLMQ KRLAMESNKI LQEQQEMERQ RRKEIAQKAA EENERYRKEM EQIVEEEEFK KKQWEEDWGS KEQLLLPKTI TAEVHPVPLR KPKYDQGVPE EEPADDLDG GTEEQGEQDF RKYEEGFDPY SMFTPEQIMG KDVRLRIKK EGSLDLAEG GVDSPIGKVV VSAVYERGAA ERHGGIVKGD EIMAINGKIV TDYTLAEAEA ALQKAWNQGG DWIDLVVAVC PPKEYDDELT

image



Background

Alternative Names	Harmonin (Antigen NY-CO-38/NY-CO-37) (Autoimmune enteropathy-related antigen AIE-75) (Protein PDZ-73) (Renal carcinoma antigen NY-REN-3) (Usher syndrome type-1C protein)
--------------------------	---

Product Name: Recombinant Human USH1C Protein
Catalog #: PEH8054



Note

For Research Use Only , Not for Diagnostic Use