



Product Name: Recombinant Human WASF2 Protein
Catalog #: PEH8055



Summary

Name	WASF2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y6W5
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FASRVSSLAE RVDRLQVKVT QLDPKEEEVS LQGINTRKAF RSSTIQDQKL FDRNSLPVPV LETYNTCDTP PPLNNLTPIYR DDGKEALKFY TDPSYFFDLW KEKMLQDTKD IMKEKRKHRK EKKDNPNRGN VNPRKIKTRK EEWEKMKMGQ EFVESKEKLG TSGYPPTLVY QNGSIGCVEN VDASSYPPPP QSDSASSPSP SFSEDNLPPP PAEFSYPVDN QRGSGLAGPK RSSVVSPSHP PPAPPLGSPP GPKPGFAPPP

image



Background

Alternative Names	Wiskott-Aldrich syndrome protein family member 2 (WASP family protein member 2) (Protein WAVE-2) (Verprolin homology domain-containing protein 2)
--------------------------	---

Product Name: Recombinant Human WASF2 Protein
Catalog #: PEH8055



Note

For Research Use Only , Not for Diagnostic Use