


Product Name: Recombinant Human PLCD1 Protein
Catalog #: PEH8056



Summary

Name	PLCD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P51178
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HGLQDDEDLQ ALLKGSQLLK VKSSSWRRER FYKLQEDCKT IWQESRKVMR TPESQLFSIE DIQEVRMGHR TEGLEKFARD VPEDRCFSIV FKDQRNTLDL IAPSPADAQH WVLGLHKIIH HSGSMDQRQK LQHWIHSCLR KADKNKDNKM SFKELQNFLK ELNIQVDDSY ARKIFRECDH SQTDSLEDEE IEAFYKMLTQ RVEIDRTFAE

image



Background

Alternative Names	1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase delta-1 (EC 3.1.4.11) (Phosphoinositide phospholipase C-delta-1) (Phospholipase C-III) (PLC-III) (Phospholipase C-delta-1) (PLC-delta-1)
--------------------------	--

Product Name: Recombinant Human PLCD1 Protein
Catalog #: PEH8056



Note

For Research Use Only , Not for Diagnostic Use