


Product Name: Recombinant Human ALDH8A1 Protein
Catalog #: PEH8057



Summary

Name	AL8A1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H2A2
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	STGEVYCRVP NSGKDEIEAA VKAAREAFPS WSSRSPQERS RVLNQVADLL EQSLEEFQA ESKDQGKTLA

image



Background

Alternative Names	2-aminomuconic semialdehyde dehydrogenase (EC 1.2.1.32) (Aldehyde dehydrogenase 12) (Aldehyde dehydrogenase family 8 member A1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use