



**Product Name: Recombinant Human CNP Protein**  
**Catalog #: PEH8060**



## Summary

<b>Name</b>	CN37_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P09543
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	46 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	MSSSGAKDKP ELQFPFLQDE DTVATLLECK TLFILRGLPG SGKSTLARVI VDKYRDGTKM VSADAYKITP GARGAFSEY KRLDEDLAAY CRRDIRILV LDDTNHERER LEQLFEMADQ YQYQVVLVPE KTAWRLDCAQ LKEKNQWQLS ADDLKCLKPG LEKDFLPLYF GWFLTCKKSE TLRKAGQVFL EELGNHKAFK KELRQFVPGD EPREKMDLVT YFGKRPPGVL HCTTKFCDYG KAPGAEYEQ QDVLKKSYSK AFTLTISALF VTPKTTGARV ELSEQQQLLW PSDVDKLSPT DNLPRGSR AH ITLGCAADVE AVQTGLDLLE ILRQEKGGSR GEEVGELSRG KLYSLGNGRW MLTLAKNMEV RAIFTGYYGK GKPVPQTQGSR KGGALQSCTI

## image



## Background

**Alternative Names** 2',3'-cyclic-nucleotide 3'-phosphodiesterase (CNP) (CNPase) (EC 3.1.4.37)

**Product Name: Recombinant Human CNP Protein**  
**Catalog #: PEH8060**

---



**Note**

For Research Use Only , Not for Diagnostic Use