


Summary

Name	HIBCH_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6NVY1
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KLAMLEEDLL ALKSPSKENI ASVLENYHTE SKIDRDKFSI LEEHMDKINS CFSANTVEEI IENLQQDGSS FALEQLKVIN KMSPTSLKIT LRQLMEGSSK TLQEVLTMEY RLSQACMRGH DFHEGVRAVL IDKDQSPKWK PADLKEVTEE DLNNHFKSLG

image



Background

Alternative Names	3-hydroxyisobutyryl-CoA hydrolase, mitochondrial (EC 3.1.2.4) (3-hydroxyisobutyryl-coenzyme A hydrolase) (HIB-CoA hydrolase) (HIBYL-CoA-H)
--------------------------	--

Note

Product Name: Recombinant Human HIBCH Protein
Catalog #: PEH8062



For Research Use Only , Not for Diagnostic Use