



Product Name: Recombinant Human PRKAG2 Protein
Catalog #: PEH8065



Summary

Name	AAKG2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UGJ0
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ATSPGGIRFF SRSRKTSGLS SSPSTPTQVT KQHTFPLESY KHEPERLENR IYASSPPDT GQRFCPSSFQ SPTRPPLASP

image



Background

Alternative Names 5'-AMP-activated protein kinase subunit gamma-2 (AMPK gamma2) (AMPK subunit gamma-2) (H91620p)

Note

For Research Use Only , Not for Diagnostic Use