


Product Name: Recombinant Human ANP32A Protein
Catalog #: PEH8067



Summary

Name	AN32A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P39687
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LRNRTPSDVK ELVLDNSRSN EGKLEGLTDE FEELEFLSTI NVGLTSIANL PKLNKLEKLE LSDNRVSGGL EVLAEKCPNL THLNLSGSKI KDLSTIEPLK KLENLKSIDL FNCEVTNLND YRENVFKLLP QLTYLDGYDR DDKEAPDSDA

image



Background

Alternative Names	Acidic leucine-rich nuclear phosphoprotein 32 family member A (Acidic nuclear phosphoprotein pp32) (pp32) (Leucine-rich acidic nuclear protein) (LANP) (Mapmodulin) (Potent heat-stable protein phosphatase 2A inhibitor I1PP2A) (Putative HLA-DR-associated protein I) (PHAPI)
--------------------------	---

Note

Product Name: Recombinant Human ANP32A Protein
Catalog #: PEH8067



For Research Use Only , Not for Diagnostic Use