


Product Name: Recombinant Human ACTN2 Protein
Catalog #: PEH8071



Summary

Name	ACTN2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P35609
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HLAEKFRQKA STHETWAYGK EQILLQKDYE SASLTEVRAL LRKHEAFESD LAAHQDRVEQ IAAIAQELNE LDYHDAVNVN DRCQKICDQW DRLGTLTQKR REALERMEKL LETIDQLHLE FAKRAAPFNN WMEGAMEDLQ DMFIVHSIEE IQLLITAHEQ FKATLPEADG ERQSIMAIQN EVEKVIQSYN IRISSNPYS TVTMDLRLTK WDKVKQLVPI RDQSLQEELA

image



Background

Alternative Names Alpha-actinin-2 (Alpha-actinin skeletal muscle isoform 2) (F-actin cross-linking protein)

Note

Product Name: Recombinant Human ACTN2 Protein
Catalog #: PEH8071



For Research Use Only , Not for Diagnostic Use