


Product Name: Recombinant Human AIMP2 Protein
Catalog #: PEH8074



Summary

Name	AIMP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q13155
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GGGAPLRVEL PTCMYRLPNV HGRSYGPAPG AGHVQEESEN SLQALESRQD DILKRLYELK AAVDGLSKMI QTPDADLDVT NIIQADEPTT LTTNALDLNS VLGKDYGALK DIVINANPAS PPLSLLVLHR LLCEHFRVLS TVHTHSSVKS VPENLLKCFG EQNKKQPRQD

image



Background

Alternative Names	Aminoacyl tRNA synthase complex-interacting multifunctional protein 2 (Multisynthase complex auxiliary component p38) (Protein JTV-1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human AIMP2 Protein
Catalog #: PEH8074

