


Product Name: Recombinant Human ATP5F1C Protein
Catalog #: PEH8080



Summary

Name	ATPG_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P36542
Host	E.coli
Species	Human
Molecular Mass	70 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VAAAKYARAE RELKPARIYG LGSLALYEKA DIKGPEDKKK HLLIGVSSDR GLCGAIHSSI AKQMKSEVAT LTAAGKEVML VGIGDKIRGI LYRTHSDQFL VAFKEVGRKP PTFGDASVIA LELLNSGYEF DEGSIIFNKF RSVISYKTEE KPIFSLNTVA SADSMSIYDD IDADVLQNYQ EYNLANIIYY SLKESTTSEQ SARMTAMDNA SKNASEMIDK LTLTFNRTRQ AVITKELIEI

image



Background

Alternative Names	ATP synthase subunit gamma, mitochondrial (ATP synthase F1 subunit gamma) (F-ATPase gamma subunit)
--------------------------	--

Note

Product Name: Recombinant Human ATP5F1C Protein
Catalog #: PEH8080



For Research Use Only , Not for Diagnostic Use