


Product Name: Recombinant Human ATP5F1C Protein
Catalog #: PEH8082



Summary

Name	ATPG_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	P36542
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LELLNSGYEF DEGSIIFNKF RSVISYKTEE KPIFSLNTVA SADSMSIYDD IDADVLQNYQ EYNLANIIYY SLKESTTSEQ SARMTAMDNA

image



Background

Alternative Names ATP synthase subunit gamma, mitochondrial (ATP synthase F1 subunit gamma) (F-ATPase gamma subunit)

Note

For Research Use Only , Not for Diagnostic Use