


Product Name: Recombinant Human ABCB6 Protein
Catalog #: PEH8083



Summary

Name	ABCB6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NP58
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MENMFDLLKE ETEVKDLPGA GPLRFQKGRI EFENVHFSYA DGRETLDQVS FTVMGQTLA LVGPSGAGKS TILRLFRFY DISSGCIRID GQDISQVTQA SLRSHIGVVP QDTVLFNDTI ADNIRYGRVT AGNDEVEAAA QAAGIHDAIM AFPEGYRTQV GERGLKLSGG EKQRVAIART ILKAPGIILL DEATSALDTS NERAIQASLA KVCANRTTIV VAHRLSTVVN ADQILVIKDG CIVERGRHEA LLSRGGVYAD MWQLQQGQEE TSEDTPQTM

image



Background

Alternative Names	ATP-binding cassette sub-family B member 6, mitochondrial (Mitochondrial ABC transporter 3) (Mt-ABC transporter 3) (P-glycoprotein-related protein) (Ubiquitously-expressed mammalian ABC half transporter)
--------------------------	---

Product Name: Recombinant Human ABCB6 Protein
Catalog #: PEH8083



Note

For Research Use Only , Not for Diagnostic Use