


Product Name: Recombinant Human RECQL Protein
Catalog #: PEH8084



Summary

Name	RECQ1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P46063
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EYCRDLIKIL KQAEELNEKL TPLKLIDSWM GKGAALKRVA GVVAPTLPRE DLEKIIAHFL IQQYLKEDYS FTAYATISYL KIGPKANLLN NEAHAITMQV TKSTQNSFRA ESSQTCHSEQ GDKKMEEKNS GNFQKKAANM LQQSGSKNTG

image



Background

Alternative Names	ATP-dependent DNA helicase Q1 (EC 3.6.4.12) (DNA helicase, RecQ-like type 1) (RecQ1) (DNA-dependent ATPase Q1) (RecQ protein-like 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use