


Product Name: Recombinant Human GET3 Protein
Catalog #: PEH8085



Summary

Name	GET3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43681
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VEPLEPTLSN IIEQRSLKWI FVGGKGGVGK TTCSCSLAVQ LSKGRESVLI ISTDPAHNIS DAFDQKFSKV PTKVKGYDNL FAMEIDPSLG VAELPDEFFE EDNMLSMGKK MMQEAMSAFP GIDEAMSYAE VMRLVKGMNF SVVVFDTAPT GHTLRLLNFP TIVERGLGRL MQIKNQISPF ISQMCNMLGL GDMNADQLAS KLEETLPVIR SVSEQFKDPE QTTFICVCIA EFLSLYETER LIQELAKCKI DTHNIIVNQL VFPDPEKPCK MCEARHKIQA KYLDQMEDLY EDFHIVKLPL LPHEVRGADK VNTFSALLLE

image



Background

Alternative Names	ATPase GET3 (EC 3.6.-.-) (Arsenical pump-driving ATPase) (Arsenite-stimulated ATPase) (Guided entry of tail-anchored proteins factor 3, ATPase)
--------------------------	---

Product Name: Recombinant Human GET3 Protein
Catalog #: PEH8085



(Transmembrane domain recognition complex 40 kDa ATPase subunit)
(hARSA-I) (hASNA-I)

Note

For Research Use Only , Not for Diagnostic Use