

Summary

Name	BCLF1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NYF8
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HSHSIQHSPE RSGSGSVGNG SSRYSPSQNS PIHHIPSRRS PAKTIAPQNA PRDESRGRSS FYPDGGDQET AKTGKFLKRF TDEESRVFLL DRGNTRDKEA SKEKGSEKGR AEGEWEDQEA LDYFSDKESG KQKFNDSEGD DTEETEDYRQ FRKSVLADQG KSFATASHRN TEEGLKYKS KVSLKGNRES DGFREEKNYK LKETGYVVER PSTTKDKHKE EDKNSERITV KKETQSPEQV KSEKLKDLFD YSPPLHKNLD AREKSTFREE SPLRIKMIAS DSHRPEVKLK MAPVPLDDSN RPASLTKDRL LASTLVHSVK

image



SDS



WB

Background

Alternative Names	Bcl-2-associated transcription factor 1 (Btf) (BCLAF1 and THRAP3 family member 1)
--------------------------	---

Product Name: Recombinant Human BCLAF1 Protein
Catalog #: PEH8086



Note

For Research Use Only , Not for Diagnostic Use