


Product Name: Recombinant Human B3GALT2 Protein
Catalog #: PEH8087



Summary

Name	B3GT2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43825
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TYTFRGFRST KSETNHSSLR NIWKETVPQT LRPQTATNSN NTDLSPQGVT GLENTLSANG SIYNEKGTGH PNSYHFKYII NEPEKCQEKs

image



Background

Alternative Names	Beta-1,3-galactosyltransferase 2 (Beta-1,3-GalTase 2) (Beta3Gal-T2) (Beta3GalT2) (EC 2.4.1.86) (UDP-galactose:2-acetamido-2-deoxy-D-glucose 3beta-galactosyltransferase 2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use