


Product Name: Recombinant Human SPAG9 Protein
Catalog #: PEH8091



Summary

Name	JIP4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	O60271
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RSHTSLKVSN SPEPQKAVEQ EDELSDVSG GSKATTPAST ANSDVATIPT DTPLKEENEG FVKVTDAPNK SEISKHIEVQ VAQETRVNST GSAENEEKSE

image



Background

Alternative Names	C-Jun-amino-terminal kinase-interacting protein 4 (JIP-4) (JNK-interacting protein 4) (Cancer/testis antigen 89) (CT89) (Human lung cancer oncogene 6 protein) (HLC-6) (JNK-associated leucine-zipper protein) (JLP) (Mitogen-activated protein kinase 8-interacting protein 4) (Proliferation-inducing protein 6) (Protein highly expressed in testis) (PHET) (Sperm surface protein) (Sperm-associated antigen 9) (Sperm-specific protein) (Sunday driver 1)
--------------------------	--

Product Name: Recombinant Human SPAG9 Protein
Catalog #: PEH8091



Note

For Research Use Only , Not for Diagnostic Use