


**Product Name: Recombinant Human PDE1C Protein**  
**Catalog #: PEH8092**



## Summary

<b>Name</b>	PDE1C_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q14123
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	16 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	EAKSQAEEGA SGKAEKKTSG ETKNQVNGTR ANKSDNPRGK NSKAEKSSGE QQQNGDFKDG KNKTDKKDHS NIGNDSKKTG GTKQRSHGSP APSTSSTCRL TLPVIKPLR HFKRPAYASS SYAPSVSKKT

## image



## Background

<b>Alternative Names</b>	Calcium/calmodulin-dependent 3',5'-cyclic nucleotide phosphodiesterase 1C (Cam-PDE 1C) (3',5'-cyclic-AMP phosphodiesterase) (EC 3.1.4.-) (3',5'-cyclic-GMP phosphodiesterase) (EC 3.1.4.35) (Hcam3)
--------------------------	---

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human PDE1C Protein**  
**Catalog #: PEH8092**

---

