


Product Name: Recombinant Human CLGN Protein
Catalog #: PEH8093



Summary

Name	CLGN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14967
Host	E.coli
Species	Human
Molecular Mass	51 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FEENSEEIDV NESELSSEIK YKTPQPIGEV YFAETFDSEGR LAGWVLSKAK KDDMDEEISI YDGRWEIEEL KENQVPGDRG LVLKSRKHH AISAVLAKPF IFADKPLIVQ YEVSFQDGID CGGAYIKLLA DTDDLILENF YDKTSYIIMF GPDKCGEDYK LHFIFRHKHP KTGVFEEKHA KPPDVLKKF FTDRKTHLYT LVMNPDDTFE VLVDQTVVVK GSLLDVPVP IKPPKEIEDP NDKKPEEWDE RAKIPDPSAV KPEDWDESEP AQIEDSSVVK PAGWLDDEPK FIPDPNAEKP DDWNEDTDGE WEAPQILNPA CRIGCGEWKP PMIDNPKYKG VWRPPLVDNP NYQGIWSPRK IPNPDYFEDD HPFLTSFSA LGLELWSMST DIYFDNFIIC SEKEVADHWA ADGWRWKIMI

image



Background

Product Name: Recombinant Human CLGN Protein
Catalog #: PEH8093



Alternative Names

Calmegin

Note

For Research Use Only , Not for Diagnostic Use