



## Summary

<b>Name</b>	CENPE_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q02224
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	45 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	SATISKEWEQ DLKSLKEKNE KLFKNYQTLK TSLASGAQVN PTTQDNKNPH VTSRATQLTT EKIRELENSL HEAKESAMHK ESKIIMQKQE LEVTNDIIAK LQAKVHESNK CLEKTKETIQ VLQDKVALGA KPYKEEIEDL KMKLVKIDLE KMKNAKEFEK EISATKATVE YQKEVIRLLR ENLRRSQQAQ DTSVISEHTD PQPSNKPLTC GGGSGIVQNT KALILKSEHI RLEKEISKLK QQNEQLIKQK NELLSNNQHL SNEVKTWKER TLKREAHKQV TCENSPKSPK VTGTASKKKQ ITPSQCKERN LQDPVPKESP KSCFFDSRSK SLPSPHPVRY FDNSSLGLCP EVQNAGAESV DSQPGPWHAS SGKDVPECKT

## image



## Background

**Alternative Names** Centromere-associated protein E (Centromere protein E) (CENP-E)

**Product Name: Recombinant Human CENPE Protein**  
**Catalog #: PEH8099**



---

(Kinesin-7) (Kinesin-related protein CENPE)

**Note**

For Research Use Only , Not for Diagnostic Use