


Product Name: Recombinant Human CSH2 Protein
Catalog #: PEH8101



Summary

Name	CSH2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P0DML3
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YIPKDQKYSF LHSQTSFCF SDIPTPSNM EETQQKSNLE LLRISLLIE SWLEPVRFLR SMFANNLVYD TSDSDDYHLL KDLEEGIQTL MGRLEDGSRR TGQILKQTYS KFDTNSHNHD ALLKNYGLLY CFRKDMDKVE

image



Background

Alternative Names	Chorionic somatomammotropin hormone 2 (Choriomammotropin) (Lactogen) (Placental lactogen) (PL)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use